

BEST PRACTICE GUIDELINES FOR JETPEEL® TREATMENTS DURING THE COVID-19 PANDEMIC



When performing a JetPeel treatment we advise to take care of your health and protect others by doing the following:

Confirm your customer does not have a fever, cough, cold, difficulty breathing or flu-like symptoms, and that they have not been in contact in the past 14 day with anyone suspected to have COVID-19.

If any of the above or true, we recommend postponing treatment until a later date.

PRIOR TO TREATMENT

Practice good hygiene to protect against infection and prevent the virus spreading.



Disinfect the treatment bed and other close surfaces with 70% alcohol and change the towels and bed cot cover between treatments.



Swipe the handpiece, tubes and device with 70% alcohol. Leave for 1 minute prior to use.

DURING THE TREATMENT



Personal protective equipment (PPE) is required. The use of gloves, a N95 facial mask, protective eye wear and a face shield/visor mask, as well as a surgical cap, shoe covers, and a fluid repellent gown are recommended.




Before putting on a mask, clean hands and forearms with alcohol-based hand rub or soap and water.



Wipe the clients face with 70% alcohol and wait 1-minute prior starting the treatment.



Lip treatments require that patients be asked to rinse the mouth with 1.5% hydrogen peroxide or 0.2% povidone-iodine for 1-minute right before the procedure.



Continue treatment as normal.



PPE is to be swapped or wiped clean between treatments (according to WHO guidelines). Fluid repellent face shields and eye protection can be used for a session of work rather than a single patient contact.

It is important to stay aware of the latest information on the COVID-19 outbreak, available through your national and local public health authority and closely follow their guidelines.