



Scalp and Hair Treatment

Presentation by: Laura Masini, International Trainer, TavTech

March 2022





Hair Matters

It is the first thing most people see when they look at another person. Taking care of your hair and scalp is just as important as caring for the rest of the skin on your body.



Scalp Matters

A healthy scalp is pink in tone, well hydrated, smooth, and free of flakes, redness, irritation, dryness, and oily. Scalp condition impacts hair growth.



Unhealthy Scalp Signs

Itching, burning, extremely oily, redness, flaky and hair loss



Use JetPeel to regenerate hair growth and rebalance the scalp

- JetPeel by TavTech is an innovative and unique aesthetics jet pressure energy-based device focused on anti-aging and skin rejuvenation, and is also used to treat **scalp and hair**.
- JetPeel is a needle-free, no pain, no downtime device using a jet stream of microdroplets to penetrate the scalp and regenerate hair growth.



Regenerate hair growth with JetPeel

- JetPeel allows for the treating of stressed, fine and damaged hair, as well as itchy, dry, oily and flaky scalp, and weakened hair structure.
- It is indicated for all skin types, and is suitable all year-round.
- This treatment is ideal also for the prevention of hair-loss, and to reinforce the hair structure from oxidative stress and from the effects of ageing.




Scalp and hair treatments

- Scalp & hair treatment with JetPeel cleanses and purifies the scalp and reduces hair-loss, while promoting hair-growth by nourishing the scalp.
- JetPeel Treatment is the best way to clear the scalp and reset the sebum levels.
- Gentle exfoliation of dead skin cells helps to unclog hair follicles, improve micro-circulation, stimulate the growth of new hair and strengthen existing hair.

Understanding the scalp types



- Before starting with the JetPeel treatment, remember to evaluate your patient's scalp condition during the consultation.
- Slightly rub the scalp with your fingers to identify sebum production. Also make sure to check for flakiness, redness, scaling and any other signs of scalp imbalances.



Dry scalp

- Typically report having dry skin overall.
- Show small white dead skin cells. Patients may complain about occasional itching and redness, as well as a feeling of tension.

Possible Causes:

- External influences and environmental factors such as pollution, air conditioning and heating.
- Chemical treatments or insufficient protection from UV rays.



Oily scalp

- Excess sebum; the scalp can look oily.
- Hair looks slick, heavy and stringy.

Possible Causes:

- Hypersecretion of the sebaceous glands
- Stress, imbalanced nourishment, hormonal disorder, medication, chemical treatments or environmental factors.



Flaky scalp

- White flakes of dead skin appear in the hair and on the scalp.

Possible Causes:

- External environmental factors.
- State of health, nutritional and psychogenic factors.

Sensitive & itchy scalp

Unpleasant sensations (tingling, burning, pain, pruritus)



Hair loss

- The hair is thinning, lacks strength, volume and density.

Possible Causes:


- Hormonal imbalance
- Stress
- Illness
- Age

JetPeel Scalp Treatment Sequence

STEP 1

MASSAGE & EXFOLIATION


Deep cleansing of the scalp; remove dead cells, excessive oil, sweat and grime; improve blood flow to strengthen the hair follicles.



STEP 2

INFUSION

This step provides effective and deep delivery of a valuable nourishment into the scalp which prevents hair loss, stimulates hair growth, and makes the hair look lush, glowing and bright



Recommended to be used with **NEW DOUBLE Handpiece**

The infusion step is performed with the NEW and INNOVATIVE DOUBLE handpiece has a double nozzle design, suitable for successful trans-epidermal hair and scalp treatments.

STEP 1

MASSAGE & EXFOLIATION



JetCare Hydro

Hyaluronic Acid and **Organic Aloe Vera**



5% Glycolic Acid-
Gentle Exfoliation

Glycolic Acid
Organic Aloe Vera and
Pro-Vitamin B₅ Panthenol

5% Glycolic acid – ideal for delicate and normal scalp



10% Glycolic Acid-
Medium Exfoliation

Glycolic Acid
Organic Aloe Vera and
Pro-Vitamin B₅ Panthenol

10% Glycolic acid – ideal for oily to flaky scalp

STEP 2

INFUSION



Scalp & Hair

Root Biotech™ - deeply nourishing for a fuller, healthy look

Vitamin B combination - niacinamide, panthenol, biotin

Redensyl™ - amino acids and antioxidants for healthier hair

Capixyl™ - a Biomimetic peptide with red clover extract (Biochanin A) which helps create thicker, fuller and healthier hair



Hydrate

Maximum concentration of **PrimalHyal™ Ultrafiller**, an innovative HA derivative

Niacinamide and **Betaine**



Regenerate

Vitamins A and E

Pentavitin® (natural plant extract)

Ragu-age®

Roman Chamomile Extract



Tight

B₃, B₅, B₇

Scalp & hair Treatment Guidelines



SCALP & HAIR

Recommended Treatment Cycle



10 Treatments

once every week for 10 weeks

MAINTENANCE (EVERY 3-4 WEEKS)

Scalp & Hair Characteristics

Itchy, Sensitive Scalp

The main reason for itchiness of scalp is inflammation and sensitivity. This may be due to environmental causes, such as pollution or sun exposure. Itchy scalp can also be caused by emotional issues such as stress and anxiety. In some cases, itchiness is caused by buildup of dry skin or from the use of harsh products that irritate the delicate skin of the scalp.

Dry Scalp

Dry scalp is sometimes genetically caused, and those who suffer from dry scalp are also often prone to other skin conditions, such as eczema. In other cases, dry skin of the scalp is caused by external or emotional factors, such as indoor heating or stress. The scalp has numerous oil glands, however, those with a tendency towards dry skin often suffer from dry scalp. Dry scalp may create a sensation of tightness and may also be accompanied by unwanted flaking.

Oily Scalp

Skin and scalp have a Hydrolipidic film consisting of sebum, sweat and moisture that forms a natural barrier to protect the skin from external elements. In those prone to oily skin, the Hydrolipidic film may over-perform, causing oiliness of the scalp. Oily skin can result from stress, hormonal changes, exposure to pollution or harsh products. An oily scalp can cause the hair roots to become greasy very quickly, which may cause acne spots.

Flaky Scalp

Flaky scalp is highly visible, as it produces small dry flakes that sit on the scalp or fall onto the clothing. The causes of flaky scalp are varied. It may occur in dry skin, oily and irritated skin, or sensitive skin. Flakiness may also be a result of contact dermatitis, which is caused by sensitivity to products used on the hair and scalp.



Scalp and Hair Treatment Guidelines

SCALP TYPE	1st Treatment	2nd Treatment	3rd Treatment	4th Treatment	5th Treatment
Itchy, Sensitive	Treatment Guideline 33 with JetCare Hydro with Boost - Regenerate	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 35 with JetCare Hydro with Boost - Hydrate
Dry	Treatment Guideline 33 with JetCare Hydro with Boost - Regenerate	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate
Oily	Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate	Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 40 with Glycolic Acid 10% with Boost - Hydrate
Flaky	Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate	Treatment Guideline 33 with JetCare Hydro with Boost - Regenerate	Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 40 with Glycolic Acid 10% with Boost - Hydrate

6th Treatment	7th Treatment	8th Treatment	9th Treatment	10th Treatment
Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 36 with JetCare Hydro with Boost - Tight	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair
Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 38 with Glycolic Acid 5% with Boost - Hydrate	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair
Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate	Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair
Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair	Treatment Guideline 41 with Glycolic Acid 5% with Scalp & Hair	Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair	Treatment Guideline 34 with JetCare Hydro with Scalp & Hair

Treatment Guideline No. 33

	MASSAGE AND EXFOLIATION	INFUSION
 Handpiece	Triple	Double
 Distance from skin	10 mm	5-7 mm
 Handpiece Angle	45-90°	90°
 Time	10 minutes	5-7 minutes
 Use	~12ml of JetCare Hydro	~6-7ml of Boost - Regenerate

Treatment Guideline No. 34

	MASSAGE AND EXFOLIATION	INFUSION
 Handpiece	Triple	Double
 Distance from skin	10 mm	5-7 mm
 Handpiece Angle	45-90°	90°
 Time	10 minutes	5-7 minutes
 Use	~12ml of JetCare Hydro	~6-7ml of of Selective Care Scalp & Hair

BEFORE




AFTER



Immediately after 10th treatment

AGE: **24**
Courtesy of Laura Masini, TavTech



BEFORE




AFTER



Immediately after 10th treatment

AGE: **24**

Courtesy of Laura Masini, TavTech



BEFORE

AFTER



Immediately after 10th treatment

AGE: 61
Courtesy of Dr. Alex Ginzburg, Israel



AGE: **24**
Courtesy of Laura Masini, TavTech



Immediately after 10th treatment

AGE: **24**
Courtesy of Laura Masini, TavTech





JETPEEL
by TavTech

Thank You

