

JetPeel by TavTech Academy



We provide people with an
opportunity to become the best
version of themselves

**We believe in
making people
feel good ...**



... by delivering
a powerful,
instant-result
treatment

Suitable for all skin types
and all year-round





JetPeel by TavTech Training Program

1. JetPeel Technology and Solutions: Training Introduction
2. Lymphatic Massage with JetPeel
 - 2.1. Lymphatic Massage with JetPeel – Face & Neck
 - 2.2. Lymphatic Massage with JetPeel – Body
3. Exfoliation with JetPeel
 - 3.1. Exfoliation with JetPeel – Face & Neck
4. Infusion with JetPeel
 - 4.1. Infusion with JetPeel – Face & Neck
5. Scalp & Hair with JetPeel
6. JetPeel Handpieces





5. Scalp & Hair with JetPeel



The JetPeel jet pressure energy system provides trans-epidermal infusion of a rich variety of nutrients and solutions

For the treatment of an array of skincare and aesthetics conditions




Anti-Aging Care



Scalp & Hair



Lightening



Clear



Soothing




Growth Factor



Renewal Care


Trans-epidermal infusion delivery system

Utilizing jet propulsion technology for needle-less, pain-free, non-invasive, non-ablative, trans-epidermal delivery



TripleJet

3 fine nozzles



DoubleJet

2 fine nozzles



SingleJet Narrow (Magic)

1 fine nozzle

JetCare - A variety of solutions, vitamins, minerals and acids




and an array of treatment guidelines

The enduring way to:

- Perform a Lymphatic Massage
- Exfoliate
- Infuse Skin Nutrients without Needles
- Stimulate Collagen and Elastin Fibers
- Tighten
- Brighten
- Eliminate Blemishes
- Smooth Skin Imperfections



JetPeel Handpieces for the Scalp & Hair Treatment



TripleJet

3 fine nozzles for
Massage and
Exfoliation



DoubleJet

2 fine nozzles for
Infusion

The Scalp and Hair Treatment Sequence

1




Recommended to be used with
TripleJet Handpiece

Massage and exfoliation with TripleJet handpiece

This step cleanses, purifies and gently exfoliates the scalp in preparation for the nourishment from the following infusion step


2



Recommended to be used with
DoubleJet Handpiece

Infusion with DoubleJet Handpiece

Deep, effective delivery of active ingredients into the scalp



The Scalp and Hair Treatment Sequence

1



Recommended to be used with
TripleJet Handpiece


Massage and exfoliation with TripleJet handpiece

This step cleanses, purifies and gently exfoliates the scalp in preparation for the nourishment from the following infusion step

2


Infusion with DoubleJet Handpiece

Deep, effective delivery of active ingredients into the scalp





- All JetPeel sessions begin with JetCare Hydro
- Gentle solution formulated for the initial massage, and to prep the skin for the deeper JetCare treatments
- Hyaluronic Acid concentration reformulated, significantly higher than in former version
- Enriched with higher concentrations of organic concentrated Aloe Vera
- Solid plastic bottle, won't puncture during transportation
- 300 ml - 20% more than in former version




Recommended to be used with TripleJet Handpiece




Divide hair by sections

1. Divide the scalp into two parts
2. Work first on one side, then repeat on the other side
3. Create lines in the hair by using a comb on long hair, or finger tips on short hair
4. Perform massage & exfoliation long these lines




Massage & exfoliation technique

1. Divide the scalp into two parts
2. Work first on one side, then repeat on the other side
3. Use the TripleJet Handpiece with Jetcare Hydro or Glycolic Acid according to Treatment Guidelines
4. Massage: The distance between the handpiece and the scalp should be 8-10mm and the handpiece angle should be held between 45-90° angle relative to the scalp
5. Exfoliation: The distance between the handpiece and the scalp should be 8-10mm and the handpiece angle should be held between 45° angle relative to the scalp
6. Move the handpiece in zigzag motions, covering the entire treatment area
7. Use a comb to create your lines in the hair on long hair
8. Start in the centre of the head (point 1), and slowly continue with points 2, 3, 4, 5, 6, 7, 8, 9, 10 and 11 until you arrive behind the patient's ear
9. Repeat on the opposite side of the head





- The ultimate exfoliating treatment line enriched with AHA's and BHA's
- The Renewal Care Series is adjusted to pH of 3.5, the ideal range for rejuvenating the skin, activating skin enzymes and barrier peptides, and for balancing the Natural Moisturizing Factor (NMF)
- Personalized exfoliation containing Glycolic Acid, Organic Aloe Vera and Pro-Vitamin B5 Panthenol
- Available in 3 concentrations:
 - 5% Glycolic acid – ideal for dry, oily and flaky scalp
 - 10% Glycolic acid – ideal for oily to flaky scalp
 - 15% Glycolic acid – ideal for oily to flaky scalp



Recommended to be used with
TripleJet Handpiece




1

Massage and exfoliation with TripleJet handpiece

This step cleanses, purifies and gently exfoliates the scalp in preparation for the nourishment from the following infusion step


2



Recommended to be used with DoubleJet Handpiece

Infusion with DoubleJet Handpiece

Deep, effective delivery of active ingredients into the scalp





- For men and women
- Creates the conditions for enhanced fullness and volume
- Formulated with innovative targeted ingredients
- Rich in Biotin Vitamin B7 that helps rejuvenate the scalp
- Contains:
 - Root Biotech™ - deeply nourishing for a fuller, healthy look
 - Redensyl™ - amino acids and antioxidants for healthier hair
 - Capixyl™ - a Biomimetic peptide with red clover extract (Biochanin A) which helps create thicker, fuller and healthier hair



**Recommended to be used with
DoubleJet Handpiece**



New and improved booster formulas designed to amplify the power of JetPeel infusions



**Recommended to be used with
DoubleJet Handpiece**



Jet Boost Hydrate (HA)

- State-of-the-art HA formula for the ultimate filler effect and beautifully smooth skin
- Maximum concentration of PrimalHyal™ Ultrafiller, an innovative HA derivative

Ideal for Itchy, sensitive, dry, oily and flaky scalp



New and improved booster formulas designed to amplify the power of JetPeel infusions



**Recommended to be used with
DoubleJet Handpiece**




Regenerate (Vit A + E)

- Specially formulated to regenerate the skin with the antioxidant power
- Serum booster enriched with vitamins A and E for high performance anti-aging action

Ideal for Itchy, sensitive, dry, and flaky scalp



New and improved booster formulas designed to amplify the power of JetPeel infusions



Recommended to be used with DoubleJet Handpiece



Tight (Vit B Complex)


- Power booster delivers effective tightening and hydrating action
- Complex of 3 vitamin B derivatives - B3, B5, B7 (vs. only B5 in previous version)
- Synergetic effect for a wide range of skin rejuvenation processes

Ideal for Itchy, sensitive, dry, and flaky scalp




Divide hair by sections

1. Divide the scalp into two parts
2. Work first on one side, then repeat on the other side
3. Create lines in the hair by using a comb on long hair, or finger tips on short hair
4. Perform infusion along these lines




Infusion technique



Recommended to be used with
DoubleJet Handpiece

1. Divide the scalp into two parts
2. Work first on one side, then repeat on the other side
3. Use the DoubleJet Handpiece with the the JetCareScalp and Hair Solution or the JetCare Boost solutions according to the Treatment Guidelines
4. The distance between the handpiece and the scalp should be 8-10mm and the handpiece angle should be held between 90° angle relative to the scalp
5. Move the handpiece in zigzag motions, covering the entire treatment area
6. Use a comb to create your lines in the hair on long hair
7. Start in the centre of the head (point 1), and slowly continue with points 2, 3, 4, 5, 6, 7, 8, 9, 10 and 11 until you arrive behind the patient's ear
8. Repeat on the opposite side of the head



SCALP & HAIR

Recommended Treatment Cycle

Initial treatment cycle:

Once every week for 10 weeks

Followed by:

Maintenance treatments every 3-4 weeks



Scalp and hair characteristics

Itchy, Sensitive Scalp

The main reason for itchiness of scalp is inflammation and sensitivity. This may be due to environmental causes, such as pollution or sun exposure. Itchy scalp can also be caused by emotional issues such as stress and anxiety. In some cases, itchiness is caused by build-up of dry skin or from the use of harsh products that irritate the delicate skin of the scalp.

Dry Scalp

Dry scalp is sometimes genetically caused, and those who suffer from dry scalp are also often prone to other skin conditions, such as eczema. In other cases, dry skin of the scalp is caused by external or emotional factors, such as indoor heating or stress. The scalp has numerous oil glands, however, those with a tendency towards dry skin often suffer from dry scalp. Dry scalp may create a sensation of tightness and may also be accompanied by unwanted flaking.

Oily Scalp

Skin and scalp have a Hydrolipidic film consisting of sebum, sweat and moisture that forms a natural barrier to protect the skin from external elements. In those prone to oily skin, the Hydrolipidic film may over-perform, causing oiliness of the scalp. Oily skin can result from stress, hormonal changes, exposure to pollution or harsh products. An oily scalp can cause the hair roots to become greasy very quickly, which may cause acne spots.

Flaky Scalp

Flaky scalp is highly visible, as it produces small dry flakes that sit on the scalp or fall onto the clothing. The causes of flaky scalp are varied. It may occur in dry skin, oily and irritated skin, or sensitive skin. Flakiness may also be a result of contact dermatitis, which is caused by sensitivity to products used on the hair and scalp.






Scalp and hair treatment guidelines

| SCALP TYPE | 1st Treatment | 2nd Treatment | 3rd Treatment | 4th Treatment | 5th Treatment |
|-----------------------------|---|---|--|---|--|
| Itchy, Sensitive | Treatment Guideline 33 with JetCare Hydro with Boost - Regenerate | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 35 with JetCare Hydro with Boost - Hydrate |
| Dry | Treatment Guideline 33 with JetCare Hydro with Boost - Regenerate | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate |
| Oily | Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate | Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 40 with Glycolic Acid 10% with Boost - Hydrate |
| Flaky | Treatment Guideline 37 with Glycolic Acid 5% with Boost - Hydrate | Treatment Guideline 33 with JetCare Hydro with Boost - Regenerate | Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 40 with Glycolic Acid 10% with Boost - Hydrate |






Scalp and hair treatment guidelines

| 6th Treatment | 7th Treatment | 8th Treatment | 9th Treatment | 10th Treatment | Maintenance (Every 3-4 Weeks) |
|---|---|---|--|---|---|
| Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 36 with JetCare Hydro with Boost Tight | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Alternate Treatment Guidelines 34 and 36 |
| Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 38 with Glycolic Acid 5% with Boost Hydrate | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Alternate Treatment Guidelines 33 and 34 |
| Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 37 with Glycolic Acid 5% with Boost Hydrate | Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Alternate Treatment Guidelines 34 and 40 |
| Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Treatment Guideline 41 with Glycolic Acid 5% with Scalp & Hair | Treatment Guideline 39 with Glycolic Acid 5% with Scalp & Hair | Treatment Guideline 34 with JetCare Hydro with Scalp & Hair | Alternate Treatment Guidelines 34 and 41 |




Treatment Guideline No. 33

| | MASSAGE AND EXFOLIATION | INFUSION |
|--|-------------------------|------------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~12ml of JetCare Hydro | ~6-7ml of Boost - Regenerate |






Treatment Guideline No. 34

| | MASSAGE AND EXFOLIATION | INFUSION |
|---|-------------------------|--|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~12ml of JetCare Hydro | ~6-7ml of of Selective Care Scalp & Hair |







Treatment Guideline No. 35

| | MASSAGE AND EXFOLIATION | INFUSION |
|--|-------------------------|---------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~12ml of JetCare Hydro | ~6-7ml of Boost - Hydrate |







Treatment Guideline No. 36

| | MASSAGE AND EXFOLIATION | INFUSION |
|---|-------------------------|-------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~12ml of JetCare Hydro | ~6-7ml of Boost - Tight |







Treatment Guideline No. 37

| | MASSAGE AND EXFOLIATION | INFUSION |
|--|--|---------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~6ml of Glycolic Acid 5% | ~6-7ml of Boost - Hydrate |
|  Remarks | After a few minutes wash with ~10ml of JetCare Hydro | |







Treatment Guideline No. 38

| | MASSAGE AND EXFOLIATION | INFUSION |
|---|--|-------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~6ml of Glycolic Acid 5% | ~6-7ml of Boost - Tight |
|  Remarks | After a few minutes wash with ~10ml of JetCare Hydro | |







Treatment Guideline No. 39

| | MASSAGE AND EXFOLIATION | INFUSION |
|--|--|--|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~6ml of Glycolic Acid 5% | ~6-7ml of of Selective Care Scalp & Hair |
|  Remarks | After a few minutes wash with ~10ml of JetCare Hydro | |

Treatment Guideline No. 40

| | MASSAGE AND EXFOLIATION | INFUSION |
|---|--|---------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~6ml of Glycolic Acid 10% | ~6-7ml of Boost - Hydrate |
|  Remarks | After a few minutes wash with ~10ml of JetCare Hydro | |

Treatment Guideline No. 41

| | MASSAGE AND EXFOLIATION | INFUSION |
|---|--|------------------------------|
|  Handpiece | Triple | Double |
|  Distance from skin | 10 mm | 5-7 mm |
|  Handpiece Angle | 45-90° | 90° |
|  Time | 10 minutes | 5-7 minutes |
|  Use | ~6ml of Glycolic Acid 5%, | ~6-7ml of Boost - Regenerate |
|  Remarks | After a few minutes wash with ~10ml of JetCare Hydro | |

Hair & scalp care

Before & After




Immediately after 10th treatment

BEFORE



AFTER



AGE: **24**

Courtesy of Laura Masini, TavTech

Immediately after 10th treatment

BEFORE



AFTER



AGE: **24**

Courtesy of Laura Masini, TavTech

Immediately after 10th treatment

BEFORE



AFTER



AGE: **24**

Courtesy of Laura Masini, TavTech

Immediately after 10th treatment

BEFORE



AFTER



AGE: **24**

Courtesy of Laura Masini, TavTech

Contraindications

- Open wounds
- Severe sensitivity to cold
- Skin neoplasms
- Active cancer
- Trigeminal nerve inflammation
- Active sinusitis or laryngitis
- Active herpes
- Active bronchial asthma
- Thrombosis
- The system has not been tested or evaluated on pregnant or lactating women, nor on children and infants, therefore is not intended for use on these groups of users

JetPeel by TavTech

The perfect trio combining:



The enduring
JetPro device



The original patented
JetPeel handpieces



The premium
JetCare solutions

Designed exclusively for use together to deliver outstanding results



JetPeel by TavTech Training Program

1. JetPeel Technology and Solutions: Training Introduction

2. Lymphatic Massage with JetPeel

2.1. Lymphatic Massage with JetPeel – Face & Neck

2.2. Lymphatic Massage with JetPeel – Body

3. Exfoliation with JetPeel

3.1. Exfoliation with JetPeel – Face & Neck

4. Infusion with JetPeel

4.1. Infusion with JetPeel – Face & Neck

5. Scalp & Hair with JetPeel

6. JetPeel Handpieces





THANK YOU!



JetPeel by TavTech Academy

