



# Spring Skin Care JetPeel Treatments

Presentation by: [Laura Masini](#), Global Trainer, TavTech





Spring into  
Skin Radiance

## The transition from cold to warm weather affects the skin

- Sun exposure
- Humidity
- Temperature
- Climate differences on day to day basis

can weaken your skin and its self-protection mechanisms



## In spring the skin suffers from:

- Increased sebum production
- Sweating
- Skin breakouts
- Dehydration
- Irritation

### **JetPeel treatments in this season should focus on:**

- Stimulating cell turnover
- Reinforcing the skin against UV radiation
- Balancing skin hydration





# Improving skin micro-circulation is a key to glowing skin

## **Pleasant and relaxing trans-epidermal massage**

- Increase oxygen, vitamins and nutrients into the tissue
- Remove cellular waste
- Enhance skin appearance
- Hydrated the skin
- Sooth the skin

# Renew your Skin, Get a Fresh Look

JetPeel by TavTech's non-invasive procedures utilizing jet-pressure energy and high-velocity jet streams of air and a **blend of AHA & BHA acids** cleanse and exfoliate the skin surface.

## Benefits are:

- Unclog pores
- Assist preventing acne breakouts
- Improve skin rejuvenation
- Even skin tone
- Boost circulation



# During the transition season your skin needs special attention

Increase skin hydration with **JetCare Wide Range Hyaluronic Acid**

Regenerate the skin tissue with **JetCare Growth Factor**

Improve even skin tone with **JetCare Lightening**

Tighten the skin with **JetCare Witch Hazel & Peptides**





## Strengthen and Boost the skin to:

- Smooth skin texture
- Increase antioxidants
- Increase collagen and elastin
- Brighten
- Lift

Three dark blue curved lines, resembling a stylized 'J' or a series of connected curves, are positioned on the right side of the slide. They start from the top right and curve downwards and to the left, then curve back up and to the right. The lines are set against a light blue rectangular background.

*4 JetPeel Spring  
Treatment Sessions*



# The JetPeel Treatment Sequence

STEP  
1

## LYMPHATIC MASSAGE

Relaxing massage for detoxifying the skin and stimulating micro-circulation

STEP  
2

## EXFOLIATION

Preps the skin to receive the nourishment from the following infusion treatment

STEP  
3


## INFUSION

Deep, effective delivery of active ingredients into the skin

STEP  
4

## BOOSTER INFUSION

Post-infusion step to amplify the treatment result



# Go Fresh Up - Renew & Detoxify

**Total treatment time:** 40 minutes

**First Treatment Session**

STEP  
1

## LYMPHATIC MASSAGE

Detox, stimulate micro-circulation, re-oxygenate and purify the skin tissue.

Use **JetCare Hydro**



STEP  
2

## EXFOLIATION

Stimulate proliferation of epidermal cells, gently exfoliate, smooth and restore the skin appearance.

Use **JetCare Glycolic Acid 16%** or **Mandelic & Salicylic**



STEP  
3

## INFUSION

Stimulate skin growth for refreshing looking skin.

Use: **JetCare Growth Factor**



STEP  
4

## BOOSTER INFUSION

Micro booster infusion to support the skin health and activate self-regeneration of the skin tissue.

Use **JetCare Regenerate**



# Peel2 Glow - ReHydrate & ReVive skin cells

**Total treatment time:** 40 minutes

## Second Treatment Session

### STEP 1

#### LYMPHATIC MASSAGE

Increase blood circulation, hydration, and natural oxygenation at the cellular level.

Use **JetCare Hydro**



### STEP 2

#### EXFOLIATION

Exfoliate dull skin and reveal healthier, radiant skin underneath.

Use **JetCare Mandelic & Salicylic**



### STEP 3

#### INFUSION

Moisturise, smooth, and plump the skin texture.

Use: **JetCare Wide Range Hyaluronic Acid**



### STEP 4

#### BOOSTER INFUSION

Micro infusion of an excellent booster for a beautiful glow effect.

Use **JetCare Glow**



# Spring Into Radiance - Radiance & Even Skin Tone

**Total treatment time:** 40 minutes

## Third Treatment Session

STEP  
1

### LYMPHATIC MASSAGE

Clean, hydrate, leave the skin soft and glowing.

Use **JetCare Hydro**



STEP  
2

### EXFOLIATION

Boost skin radiance, stimulate cell renewal, improve clarity and uniformity of the skin complexion.

Use **JetCare Glycolic Acid 10%** or **Mandelic & Salicylic**



STEP  
3

### INFUSION

Effective skin lightening designed to target uneven skin tone, and improve pigmentation.

Use: **JetCare Lightening**



STEP  
4

### BOOSTER INFUSION

Protects the skin barrier and smooths skin texture.

Use **JetCare Hydrate**



# Spring Holiday Preparation

**Total treatment time:** 40 minutes

**Fourth Treatment Sessions**

STEP  
1

## LYMPHATIC MASSAGE

Boost the skin barrier, hydrate and smooth the skin.

Use **JetCare Hydro**



STEP  
2

## EXFOLIATION

Accelerate the cell turnover decongesting pores, reducing irritation and improving skin tone.

Use **JetCare Mandelic & Salicylic**



STEP  
3

## INFUSION

Reinforce and protect the skin, reduce puffiness, restore firmness for youthful and fresh looking skin.

Use: **JetCare Witch Hazel & Peptides**



STEP  
4

## BOOSTER INFUSION

Prevent cell damage and water loss, keeping the skin moist and plump.

Use **JetCare Tight**





# JetPeel by TavTech – the perfect trio combining:

- the enduring JetPro device,
- the original patented JetPeel handpieces, and
- premium JetCare solutions

**designed exclusively for use together to deliver outstanding results**





JETPEEL  
by TavTech

THANK YOU!

