



JetPeel's JetPro by TavTech Devices

The Gentle Revolution in Aesthetic Technology

Revolutionize aesthetic treatment with JetPeel by TavTech – exclusive jet pressure-energy technology that takes skincare results to the next level of beauty. JetPeel innovation is based on the JetPro by TavTech devices that generate fine jet streams of pressurized air, delivered to the skin with JetPeel's exclusive handpieces, bringing the active ingredients of JetCare by TavTech liquid solutions to the deeper skin layers for the most effective and visible impact.

JetPeel treatment is completely non-invasive, needle-free, and non-ablative. There is no skin breakage, no scarring and no recovery or downtime. Just a comfortable, pleasant experience, for the treatment of a wide range of skin concerns, for all skin types, all year round.

Choose from the range of JetPro by TavTech devices to enhance your clinic offerings and bring the JetPeel revolution to your customers:



JetPro by TavTech

The iconic JetPro is the original JetPeel device, creating fine misty jet streams that impact the skin at 95 p.s.i. and at a speed of 200 m/sec. Combined with the JetPeel handpiece and JetCare liquid solutions collection, JetPro is the ideal option for classic JetPeel treatment in the clinic or beauty salon.



JetPro Duo by TavTech

JetPro Duo combines the JetPro jet pressure technology with LED light therapy for an integrated, simultaneous Jet-LED treatment. The device features red LED add-on for an outstanding anti-aging and smoothing effect, and blue LED add-on for effective treatment of acne prone skin. Based on NASA developed technology, the JetPro Duo is the ideal option to provide enhanced JetPeel results, with no extra treatment time required.



JetPro ToGo by TavTech

Take your clinic anywhere and expand your business potential. JetPro ToGo is the portable version of JetPro Duo, a smaller, compact trolley-like device that is easily transported and is mobile. Ideal for providing dual Jet-LED treatment for bridal parties, at customers' homes, at business events and conferences, or wherever your clientele request it. Same technology, same results, and now in portable mode.

TavTech: Beauty has no borders

TavTech is the developer of JetPeel, JetCare and JetPro, the world's leading aesthetic treatment device based on jet pressure technology. TavTech products can be found in clinics globally, used by cosmetic specialists, aestheticians and dermatologists who trust in JetPeel's exceptional and visible results. Based on over two decades of R&D and a commitment to quality, TavTech continues to push the boundaries of innovation and expand its portfolio of excellence in aesthetic solutions.

